

Item and Specification:

Standard Adhesive - 100 micron tissue tape

<i>With 100 mic. Tissue Tape (Standard)</i>		<i>With 205 mic. Polyester Tape (Premium)</i>	
1.1	12.7mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNS), No Groove, MUV	2.1	12.7mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNS), No Groove, MUV
1.2	20mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNSNS), No Groove, MUV	2.2	20mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNSNS), No Groove, MUV
1.3	25.4mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNSNSNS), No Groove, MUV	2.3	25.4mm x 1.5mm x 30m 3.175mm pitch on one side (SNSNSNSNS), No Groove, MUV

Size and Tolerance:

Standard Adhesive - 100 micron tissue tape

<i>Item</i>	<i>Lenght (m)</i>	<i>Width (mm)</i>			<i>Thickness</i>	<i>Thickness of Tissue Tape</i>
Size	30	12.7	20	25.4	1.5mm	0.1
Tolerance	0/+1	±0.2	±0.2	±0.2	±0.15	±0.01

Premium Adhesive - 250 micron tissue tape

<i>Item</i>	<i>Lenght (m)</i>	<i>Width (mm)</i>			<i>Thickness</i>	<i>Thickness of Tissue Tape</i>
Size	30	12.7	20	25.4	1.5mm	0.205
Tolerance	0/+1	±0.2	±0.2	±0.2	±0.15	±0.02

Magnetic Data:

<i>Grade Code</i>	<i>BR (GS)</i>	<i>HCB (Oe)</i>	<i>HCI(Oe)</i>	<i>BH(Max)(MGOe)</i>
NF06 Extrusion	≥1600	≥1200	≥1800	0.60 - 0.80

<i>Grade Code</i>	<i>Thickness (mm)</i>	<i>Pole Pitch (mm)</i>	<i>Pull (g/cm₂)</i>
NF06 Extrusion	1.5	3.175	≥55